Fixed Diluter Air Entrainment Mask

Instructions for use

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en

Intended Use:

The Fixed Diluter Air Entrainment Mask is used to deliver fixed concentrations of oxygen to a patient during oxygen therapy.

Directions for Use:

NOTE: These directions are general guidelines intended for use by qualified medical personnel.

1. Attach the corrugated tubing to the mask grummet.
2. Select the appropriate oxygen diluter (see table).
3. Attach the diluter to the corrugated tubing.
 NOTE: Oxygen concentration and recommended flow rate are imprinted on each diluter.
4. If desired, press the high humidity adapter over the distal end of the diluter. Connect corrugated tubing (not supplied) between the adapter and a humidity source.
 WARNING: Use only room air with the humidity source. Use of oxygen will affect the desired concentration.
5. Connect supply tubing between the diluter and an appropriate oxygen source.
6. Adjust the oxygen flow to the appropriate level (see table) and check for gas flow through the device. Listed flow rates are recommendations only. Individual circumstances may require flows other than those suggested.
7. Place the mask on the patient's face with the elastic strap below the ears and around the neck. Gently pull the strap ends until the mask is secure. Mold the metal strip on the mask to fit the face.

Colour	FiO2	Recommended O2 Litre Flow
Blue	24%	4LPM
White	28%	4LPM
Orange	31%	6LPM
Yellow	35%	8LPM
Red	40%	10LPM
Pink	50%	10LPM
Green	60%	12LPM

WARNING:

1. Use of a humidifier with a low concentration diluter (i.e. 24%, 28%, 31%) may cause the humidifier pressure relief valve to activate and consequently affect gas flow to the patient.
2. Prior to use, ensure all connections are secure. Do not over tighten or force connections. Do not alter device. Improper connections may result in disconnection or leakage during use, which may result in performance degradation and increase the risk of respiratory distress
3. Do not use in the presence of an ignition source.
4. Prior to use, ensure compatibility with gas outlet. Standard connectors are compatible with 5-7mm male fittings.
5. Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
6. Do not use if package is open or damaged.
7. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Disposal Instructions:

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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