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REF 41007S
REF 1007SG
REF 1008SG

RL24SWY 2024-12

REF	Tubing Included	Length	Tubing Color	Patient Size
41007S	✓	2.1 m	Clear	Adult
1007SG	✓	2.1 m	Green	Adult
1008SG	✓	2.1 m	Green	Pediatric

Tubing Included; Schlauch enthalten; Tubulure incluse; Inclusief slang; tubo incluso; Tubos incluídos; Tubo incluído; Vrátane hadičiek; Včetně kroužku s háčkem; Slang medfølger; Putket mukana; Slangen medfølger; Slange følger med; Περιλαμβάνεται σωλήνας Rurka w zestawie; Priložena cevka; Tuburi incluse; Включени тръби; Cső mellékelve; Hortum dahildir.

Length; Lange; Longueur; Lengte; Lunghezza; Longitud; Comprimento; Dužina; Dĺžka; Delka; Langd; Pituus; L.ngde; Lengde; Μήκος; Długość; Dolžina; Lungime; Дължина; Uzunluk.

Tubing Colour – Clear, Green; Schlauch Farbe - Transparent,Grün; Couleur de tubulure - Transparent, vert, Slangkleur - doorzichtig, groen, Colore del tubo - verde chiaro, Color del tubo - Verde transparente, Color do tubagem - Transparente, Verde, Farba hadičiek – číra, zelená, Barva hadičky – čirá, zelená, Slangens färg – Klar, Grön, Letkun väri – Kirkas, Vihreä, Slangefarve – klar, grøn, Farge på slange – Klar, grønn, Χρώμα σωλήνων – Ανοιχτό, Πράσινο, Kolor rurki – przezroczysty, zielony, Barva cevke - prozorna, zelena, Culoare tubulatură – Transparent, verde, Цвят на тръбата – Прозрачен, Зелен, Tömlő színe – Átlátszó, Zöld, Tüp Rengi – Şeffaf, Yeşil.

Patient Size - Adult, Pediatric; Patienten Größe - Erwachsene, Kinder; Taille patient - Adulte, Pédiatrique; Patiëntenmaat - Volwassenen, Pediatrisch; Misura adulto, pediatrico; Tamanho - Adulto, Pediátrico; Tamanho Paciente - Adulto, Pediátrico; Veľkosť pacienta - Dospelý, Pediatrická; Velikost pacienta - Dospělý, Pediatr; Patientstorlek - Vuxen, Pediatrisk; Potilaan koko - aikuiset, lapset; Pasientstørrelse – voksen, barn; Pasientstørrelse - Voksen, Pediatrisk; Μέγεθος ασθενούς - Ενηλίκων, Παιδιατρικό; Rozmiar pacjenta – dorosły, dziecko; Velikost bolnika – odrasli, pediatrični bolniki; Dimensiune pacient - adult, pediatric; Пациентски размер - възрастни, педиатрични пациенти; Betegméret - felnőtt, gyermek; Hasta boyutu - Yetişkin, Pediatrik.

en

High Concentration Oxygen Mask with 2.1m Star Lumen Tubing and Standard Connector

Intended use:
For administration of oxygen by inhalation for patients experiencing hypoxia.

These devices are intended for use by medical personnel, lay persons, or patients that have been prescribed oxygen.

Intended to be used for a maximum of 30 days.

Instructions for use:
1. Attach supply end of tubing by pushing end connector directly onto compatible gas outlet.
2. Set oxygen to prescribed flow as directed by trained clinician.
3. Check for secure connections by confirming flow of gas through device.
4. Place the mask on the patient's face with the elastic strap below the ears and around the neck.
5. Gently pull the strap ends until the mask is secure.
6. Mold the metal strip on the mask to fit the face.

Disposal Instructions:
Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Warning:
Prior to use, ensure all connections are secure. Do not over tighten or force connections. Do not alter device. Improper connections may result in disconnection or leakage during use, which may result in performance degradation and increase the risk of respiratory distress.

Do not use in the presence of an ignition source.

Prior to use, ensure compatibility with gas outlet. Standard connectors are compatible with 5-7mm male fittings.

Do not use if device is soiled or damaged.

Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

Intended to be used for a maximum of 30 days.

Max Pressure: 50 PSI. Max Flow: 15LPM

Non-Sterile. Do not sterilize prior to use.

If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Caution should be exercised when delivering oxygen via a face mask to patients with chronic obstructive pulmonary disease or facial trauma.

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DEHP free, DEHP-frei, Sans DEHP, DEHP vrij, Senza DEHP, Sin DEHP, Sem DEHP, Bez DEHP, Neobsahuje DEHP, DEHP-fri, DEHP-vapaa, DEHP-fri, DEHP-fri, Χωρίς DEHP, Nie zawiera DEHP, Brez DEHP, Färä DEHP, не съдържа DEHP, DEHP mentes, DEHP içermez.

Ningbo Shengyurui Medical Appliances Co., Ltd.
No. 138, Binhaisi Road, Hangzhou Bay New Zone. 315336, Ningbo, Zhejiang Province, P.R CHINA

CE 1639

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

CH REP MedNet Swiss GmbH
D4 Platz 4, 6039 Root D4, Switzerland

CH importer: Medline International Switzerand Sàrl, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

UK representative: Etheria Medical Ltd, The Old brush Factory Unit 2D Whickham Industrial Estate Swalwell Newcastle Upon Tyne NE16 3DA United Kingdom
Importer in UK: Medline Industries Ltd, WA4 6HL, UK

Medline International B.V.,
Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

Distributed in Australia by:
Medline International Two Australia Pty. Ltd., Sydney, Australia