# UROSID® IAP Manometer

**IU-UROSID/00 (08/11/2023)**

CE
0123
Dir. MDD 93/42/EEC

Bioengineering Laboratories Srl
Via Provinciale per Dubino 10,
23014, Delebio (SO), Italy

CH REP

CMC Medical Devices GmbH
Bahnhofstrasse 32, CH-6300 Zug. Switzerland

CH importer: Medline International
Switzerland SARL, c/o MN & Associés SA, 1
Place de Longemalle, 1204 Genève,
Switzerland

ISTRUZIONI PER L'USO/ИНСТРУКЦИИ ЗА УПОТРЕБА/NÁVOD K POUŽITÍ/BRUGSVEJLEDNING/ANWENDUNGSHINWEISE/ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ/INSTRUCTIONS FOR USE/INSTRUCCIONES DE USO/KASUTUSJUHISED/KÄYTTÖOHJEET/INSTRUCTIONS D'EMPLOI/UPUTE ZA UPORABU/HASZNÁLATI UTASÍTÁSOK/NAUDOJIMO INSTRUKCIJOS/LIETOŠANAS INSTRUKCIJA/BRUKSINSTRUKSJONER/INSTRUCTIES VOOR GEBRUIK/INSTRUKCJE UŻYWANIA/INSTRUÇÕES DE USO/INSTRUCȚIUNI DE UTILIZARE/POKYNY PRE POUŽITIE/NAVODILA ZA UPORABO/BRUKSANVISNING/KULLANMA TALİMATLARI

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**EN - INTRA-ABDOMINAL PRESSURE MEASUREMENT SYSTEM**

**Intended use:** Foley Manometer is a sterile, disposable device for monitoring intra-abdominal pressure ('IAP') through the bladder (bladder pressure or 'P-Vesical') in hospitalised patients. (Accuracy: ± 5% of pressure reading). The Foley Manometer medical device is also intended for use in combination with Foley Catheters. While using the Foley Manometer device, it is also possible to sample urine

**Target population.**

Patients with a body weight of 10 kg or more, with a Foley catheter, and at risk of developing intra-abdominal hypertension.

**Intended users**

The Foley Manometer device must be used by personnel trained and qualified in intra-abdominal pressure ('IAP') monitoring procedures.

**CONTRAINDICATIONS**

No known contraindications for use as per INTENDED USE/ INDICATIONS FOR USE.

**INSTRUCTIONS FOR USE (refer to the drawings printed on the pouch).**

**1)Initial set-up and preparations.**

A) Follow strictly aseptic techniques.

B) Open the pouch containing the device and close the red tube clamp.

C) Place the urine collection device **UNDER** the patient's bladder and tape the drainage tube to the bed sheet to secure it in place. Thoroughly disinfect the catheter connection before disconnecting and inserting the Foley manometer device between the catheter and the drainage device.

D) Prepare the Foley Manometer device with 20ml sterile saline solution through its **needleless** sample port. **Prepare only** in the initial set-up phase or later to remove the air from the **Foley Manometer** tube.

**2)Urinary drainage:**

**1)**Leave the **Foley manometer** in its horizontal drainage position with the clamp closed, as shown in the illustrations.

**N.B.: Thoroughly disinfect the needleless sample port before taking the urine sample.**

**3)Bladder pressure monitoring:**

1) Place the '0 mmHg' mark of the tube on the mid-axillary line of the iliac crest and raise the filter vertically above the patient. Open the red clamp and read **P-Vesical** (end-expiration value) when the meniscus has stabilised.

2) Close the red clamp and place the **Foley Manometer device** in its horizontal drainage position.

**Warning**: for **anuric** patients, a new preparation is **not** necessary each time the pressure is determined. The fluid returned to the bladder from the tube provides adequate volume for reliable pressure determination.

**4) Urine sampling:**

The final urine produced by the patient is in the manometer tube, next to the Foley catheter. A urine sample can be taken

at any time by temporarily opening the red clamp and withdrawing the necessary volume from the sample port, using a 10 ml or 20 ml syringe. After sampling: close the red clamp and replace the volume taken with sterile saline solution to fill the manometer tube.

**WARNINGS:**

- Read the instructions before use.
- Open the primary packaging following the directions (PEEL).
- Intended for single-patient use only.
- Do not re-sterilise.
- DO NOT RE-USE. Reuse may cause infection or other illnesses/injuries and may impair the functionality of the device.
- Use protective gloves.
- The product must be used immediately after opening the packaging.
- In case of improper transport and/or handling, the device or packaging may be structurally and/or functionally damaged.
- Do not expose to heat or direct sunlight.
- After use, dispose of the device according to the facility's disposal policy.
- Never empty the Foley manometer tube into the urine drainage device during use.
- Always check that the Foley manometer tube is constantly filled with fluid when draining.
- A slow descent (>20 - 30 seconds) of the meniscus in the tube during a P-Vesical determination means that the Foley catheter is blocked or kinked.
- Avoid U-bends of the main urinary drainage tube (which would obstruct urinary drainage).
- Replace the Foley manometer every 7 days or sooner if the Foley catheter or urine collection device is replaced before the expiration of the 7-day maximum usage period.
- For **anuric patients**, a new preparation is **not** necessary each time the pressure is determined. The fluid returned to the bladder from the tube provides adequate volume for reliable pressure determination

**Note: Any serious incident occurring with the device must be reported to the manufacturer and to the Competent Authority of the Member State in which the user and/or patient is established.**

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