# SureSite® SELECT Dressing

REF **VAD1689-INT**

EN SureSite® SELECT Dressing
BG Обяснение на символите на опаковката
CS Legenda značek na obalu
DA SureSite® SELECT-forbinding
DE SureSite® SELECT-Verband
EL Επίδεσμος SureSite® SELECT
ES Apósito SureSite® SELECT
ET Lehel olevate sümbolite selgitus
FI SureSite® SELECT -sidos
FR Pansement SureSite® SELECT
GA Míniú na Siombailí ar an bPacáiste
HR Objašnjenje simbola na pakiranju
HU A csomagoláson található szimbólumok magyarázata
IT Cerotto per medicazione SureSite® SELECT
LT Ant pakuotės esančių simbolių paaiškinimas
LV Uz iepakojuma esošo simbolu skaidrojums
MT SureSite® SELECT Dressing
NL SureSite® SELECT-verband
PL Opatrunek SureSite® SELECT
PT Penso SureSite® SELECT
RO Pansament SureSite® SELECT
SK Obväz SureSite® SELECT
SL Obveza SureSite® SELECT
SV SureSite® SELECT förband
AR SureSite® SELECT ضمادة
JA SureSite® SELECT ドレッシング材
NB SureSite® SELECT kompress
TR SureSite® SELECT Pansuman Örtüsü

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EN **Indications:** As a primary dressing, to secure vascular access line sites. Can be left in place, if clean and dry, for up to 7 days. Keep dry during use. Remove/change at first sign of soiling or moisture.
If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.
**General Precautions & Observations:** Can be used as a primary or secondary dressing. Should any signs of redness or irritation appear, discontinue use and consult a healthcare professional.
Used devices should be disposed of in proper waste containers per hospital protocol.
Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
Distributed and imported by

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Australian Sponsor:
Emergo Australia Level 20 Tower II,
Darling Park 201 Sussex Street
Sydney, NSW 2000 Australia

MD

UK CA 0086 **UK approved body:** BSI Assurance UK Ltd, Kitemark Court Davy Avenue, Knowlhill, Milton Keynes MK5 8PP, UK (Approved body number 0086)
**UK representative and importer in UK:** Medline Industries Ltd, WA4 6HL, UK

Manufactured by;
Centurion Medical Products, LP
100 Centurion Way
Williamston, MI 48895 USA

Manufactured for Medline Industries, LP
Medline and SureSite are registered trademarks of Medline Industries, LP
Three Lakes Dr. Northfield, IL 60093 USA

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EC REP Medline International France SAS