SureSite® SELECT WR Dressing

REF **VAD1680WR-INT**

EN SureSite® SELECT WR Dressing
BG Превръзка SureSite® SELECT WR
CS Kryti SureSite® SELECT WR
DA SureSite® SELECT WR-forbinding
DE SureSite® SELECT WR-Verband
EL Επίδεσμος SureSite® SELECT WR
ES Apósito SureSite® SELECT WR
ET SureSite® SELECT WR-i side
FI SureSite® SELECT WR -sidos
FR Pansement SureSite® SELECT WR
GA Cóiriú SureSite® SELECT WR
HR Povoj SureSite® SELECT WR
HU SureSite® SELECT WR kötszer
IT Cerotto per medicazione SureSite® SELECT WR
LT „SureSite® SELECT WR“ tvarstis
LV SureSite® SELECT WR pārsējs
MT SureSite® SELECT WR Dressing
NL SureSite® SELECT WR-verband
PL Opatrunek SureSite® SELECT WR
PT Penso SureSite® SELECT WR
RO Pansament SureSite® SELECT WR
SK Obväz SureSite® SELECT WR
SL Obveza SureSite® SELECT WR
SV SureSite® SELECT WR förband
AR SureSite® SELECT WR ضمادة
JA SureSite® SELECT WR ドレッシング材
NB SureSite® SELECT WR kompress
TR SureSite® SELECT WR Pansuman Örtüsü

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EN **Indications:** As a primary dressing, to secure vascular access line sites. Can be left in place, if clean and dry, for up to 7 days. Keep dry during use. Remove/change at first sign of soiling or moisture.
If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.
General Precautions & Observations: Can be used as a primary or secondary dressing. Should any signs of redness or irritation appear, discontinue use and consult a healthcare professional.
Used devices should be disposed of in proper waste containers per hospital protocol.
Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
Distributed and imported by

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Australian Sponsor:
Emergo Australia Level 20 Tower II,
Darling Park 201 Sussex Street
Sydney, NSW 2000 Australia

MD

UK CA 0086

UK approved body: BSI Assurance UK Ltd, Kitemark Court Davy Avenue, Knowlhill, Milton Keynes MK5 8PP, UK (Approved body number 0086)
UK representative and importer in UK: Medline Industries Ltd, WA4 6HL, UK

Manufactured by;
Centurion Medical Products, LP
100 Centurion Way
Williamston, MI 48895 USA

Manufactured for Medline Industries, LP
Medline and SureSite are registered trademarks of Medline Industries, LP
Three Lakes Dr. Northfield, IL 60093 USA

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EC REP Medline International France SAS