HUDSON RCI® REF 3230

# Disposable Humidifier
## with 4 PSI Pressure Relief Valve

**Rx ONLY**

Not made with natural rubber latex.
Nicht aus Naturkautschuklatex hergestellt.
Non fabriqué en caoutchouc de latex naturel.
Non fabbricato in lattice di gomma naturale.
No fabricado con látex de caucho natural.
Não fabricado com látex de borracha natural.
Niet gemaakt met natuurlijke rubberlatex.
Inte tillverkad av naturligt latexgummi.
Ei valmistettu luonnonkumilateksista.
Ikke lavet af naturligt gummilatex.
Inneholder ikke naturgummilateks.
Δεν είναι φτιαγμένο από φυσικό ελαστικό λάτεξ.
Wyprodukowano bez użycia naturalnego lateksu.
Doğal kauçuk lateksten üretilmemiştir.

Teleflex MEDICAL
3015 Carrington Mill Blvd
Morrisville, NC 27560 USA
USA: (866) 246 6990
International: +1 919 544 8000

CE 1639
EC REP
Teleflex Medical
IDA Business and Technology Park
Dublin Rd,
Athlone, Co. Westmeath
Ireland

© Copyright 2007, 2010, 2013, 2016, 2018 Teleflex

L02855 R01

Hudson RCI is a registered trademark of Teleflex Incorporated or its affiliates.

en

## Disposable Humidifier with 4 PSI Pressure Relief Valve

**Directions for Use:**
NOTE: These directions are general guidelines intended for use by qualified medical personnel.

1. Fill the humidifier jar to the maximum fill line with sterile, distilled or bottled water. Tightly secure the lid to the jar.
2. Connect the humidifier lid wing nut to the oxygen flowmeter outlet.
3. Attach oxygen tubing (not supplied) to the humidifier lid outlet stem. Be sure all connections are secure.
4. Turn the oxygen source to the prescribed flow.

**CAUTION:** An alarm whistle from the pressure relief valve will activate if the jar accumulates too much pressure. Inspect the system for excessive blockage or a kinked tube.
To test the pressure relief valve, occlude oxygen tubing at 3 LPM. Listen for an alarm condition.

**WARNING:** Oxygen supplies of 10 psi or less using low flow (less than 2 LPM) may affect the audible alarm.

**WARNING:** Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

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