MICRO MIST® and UP-DRAFT II® Small Volume Nebulizers

Instructions for use

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Table 1: Device Information

Nebulizer Category	Nebulizer SKU	Connector	Patient Population	Flow Rates	Angles of Operation	Shared Specifications
UP-DRAFT II®	MN1705CE	Standard	Adult	5-9 LPM	0-45 degrees from vertical	Driving Gas: Compressed Air Driving Gas Pressure: 345 kPa / 50 PSI Maximum Fill Volume: 6mL Average Nebulizer Kit Mass: 0.051kg
	MN1707CE	Standard	Pediatric			
	MN1730CE	None	Adult and Pediatric			
	MN1732CE	Standard	Adult and Pediatric			
	MN1734CE	Standard	Adult and Pediatric			
MICRO MIST®	MN1880CE	None	Adult and Pediatric	8 LPM	0-90 degrees from vertical	
	MN1883CE	Standard	Adult and Pediatric			
	MN1890CE	None	Adult and Pediatric			
	MN1891CE	Standard	Adult and Pediatric			
	MN1892CE	Standard	Adult and Pediatric			
	MN1893CE	None	Adult			
	MN1894CE	Standard	Pediatric			
	MN1895CE	Universal	Adult and Pediatric			
	MN1898CE	Universal	Adult			

Table 1: Device Information; Tabelle 1: Geräteinformationen; Tableau 1 : Informations concernant le dispositif; Tabel 1: Toestelinformatie; Tabella 1: Informazioni sul dispositivo; Tabla 1: Información del dispositivo; Tabela 1: Informação sobre o dispositivo; Tabuľka 1: Informácie o zariadení; Tabulka 1: Informace o prostředku; Tabell 1: Information om enheten; Taulukko 1: Laitetiedot; Tabel 1: Oplysninger om enheden; Tabell 1: Informasjon om enhet; Πίνακας 1: Πληροφορίες συσκευής; Tabela 1: Informacje o urządzeniu; Tabela 1: Informacije o pripomočku; Tabel 1: Informații despre dispozitiv; Таблица 1: Информация за устройството; 1. táblázat: A készülékre vonatkozó információk; Tablo 1: Cihaz Bilgileri; الجدول 1: معلومات عن الجهاز

Nebulizer Category; Vernebler-Kategorie; Catégorie nébuliseur; Vernevelaar categorie; Categoria nebulizzatore; Categoría del nebulizador; Categoria do nebulizador; Kategória nebulizéra; Kategorie nebulizátoru; Kategori av nebulisator; Sumutinluokka; Forstøverkategori; Kategori av forstøver; Κατηγορία νεφελοποιητή; Kategoria nebulizatorów; Kategorija nebulizatorja; Categorie Nebulizator; Категория небулизатор; Porlasztó kategória; Nebulizatör Kategorisi; فئة جهاز البخار

Nebulizer SKU; Vernebler-SKU; Nébuliseur SKU; Vernevelaar SKU; Nebulizzatore SKU; Referencia del nebulizador; SKU do nebulizador; Nebulizér SKU; Nebulizátor SKU; Nebulisator SKU; Sumutin SKU; Forstøverens varenummer; Forstøver SKU; SKU νεφελοποιητή; SKU nebulizatora; SKU nebulizatorja; SKU Nebulizator; Небулизатор SKU; Porlasztó SKU; Nebulizatör SKU; رقم تعريف جهاز الاستنشاق

Connector - Standard, Universal, None; Konnektor - Standard, Universal, keine; Connecteur - Standard, Universel, Aucun; Verbindingsstuk - Standaard, Universeel, Geen; Connettore - Standard, Universale, Nessuna; Conector - Norma, Universal, Ninguno; Conetor - Standard, Universal, Nenhum; Konektor - Štandardné, Univerzálny, Žiadne; Konektor - Standardní, Univerzální, Žádné; Anslutning - Standard, Universell, Ingen; Liitin - Standardi, Yleinen, Ei mitään; Stik - Standard, Universal, Ingen; Koblingsstykke - Standard, Universell, Ingen; Συνδετήρας - Στάνταρ, Καθολικός, Κανένα; Przyłącze - Standardowe, Uniwersalne, Brak; Priključek - Standardno, Univerzalni, Nobeno; Conector - Standard, Universal, Nici unul; Конектор - Стандартен, Универсален, Няма; Csatlakozó - Szabványos, Univerzális, Nincs; Konektör - Standart, Evrensel, Yok; الموصل-الأساسي, العالمي, لا يوجد

Patient Population - Adult, Pediatric; Patientenpopulation - Erwachsene, Kinder; Population de patients - Adulte, Pédiatrie; Patiëntengroep - Volwassenen, Kinderen; Popolazione dei pazienti - Adulti, Bambini; Población de pacientes - Adulta, Pediátrica; População de pacientes - Adulto, Pediátrico; Populácia pacientov - Dospelí, Deti; Populace pacientů - Pro dospělé, Pro děti; Patientpopulation - Vuxen, Pediatrisk; Potilasryhmä - Aikuinen, Lapsi; Patientpopulation - Voksen, Pædiatrisk; Pasientpopulasjon - Voksen, Pediatrisk; Πληθυσμός ασθενών - Ενήλικα, Παιδιατρικό; Populacja pacjentów - Dorośli, Dzieci; Populacija bolnikov - Odrasli, Otroci; Grupul de pacienți - Adult, Copil; Пациентска популация - Възрастни, Педиатрични; Betegpopuláció - Felnőtt, Gyermekgyógyászati; Hasta Popülasyonu - Yetişkin, Pediatrik; عدد المرضى-للبالغين, للأطفال

Angles of Operation - 0-90° from vertical, 0-45° from vertical; Betriebswinkel - 0 - 90° von der Senkrechten, 0 - 45° von der Senkrechten; Angles d'utilisation - 0-90° par rapport à la verticale, 0-45° par rapport à la verticale; Werkingshoeken - 0-90° verticaal, 0-45° verticaal; Angoli di funzionamento - 0-90° da posizione verticale, 0-45° da posizione verticale; Ángulos de operación - 0-90° de la vertical, 0-45° de la vertical; Ângulos de funcionamento - 0-90° a partir da vertical, 0-45° a partir da vertical; Uhly prevádzky - 0 – 90° od vertikály, 0 – 45° od vertikály; Provozní úhly - 0-90° od svislé polohy, 0-45° od svislé polohy; Driftsvinklar - 0-90° från vertikal, 0-45° från vertikal; Toimintakulmat - 0-90° pystysuorasta, 0-45° pystysuorasta; Betjeningvinkler - 0-90 ° fra lodret, 0-45 ° fra lodret; Driftsvinkler - 0-90° fra vertikal, 0-45° fra vertikal; Γωνίες λειτουργίας - 0-90° από την κάθετη, 0-45° από την κάθετη; Kąty działania - 0-90° od pionu, 0-45° od pionu; Koti delovanja - 0–90° navpično, 0–45° navpično; Unghiuri de Operare - 0-90° din poziție verticală, 0-45° din poziție verticală; Ъгли на работа - 0-90° от вертикала, 0-45° от вертикала; Működési szögek - 0-90° függőlegesen, 0-45° függőlegesen; Çalışma Açıları - dikeyden 0-90°, dikeyden 0-45°; زوايا التشغيل- 0-90° من الوضع الرأسي, 0-45° من الوضع الرأسي

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en

Intended purpose:

Medline Hudson RCI Small Volume Nebulizers are intended to be used for intermittent aerosol therapy to deliver prescribed solutions for inhalation. Intended for use on pediatric or adult patients. Intended for use by medical personnel, lay persons, or patients that have been prescribed aerosol therapy. For use in acute and non-acute care settings. Intended for single patient use.

Directions for Use:

1. Place the prescribed solution into the nebulizer jar.
2. Attach the nebulizer cap to the nebulizer jar with a clockwise rotation.
3. If provided, connect components accordingly: attach mask to cap; attach tee to cap; attach reservoir tube to tee; attach mouthpiece to tee. Discard any items if not used.
4. Attach tubing by pushing end connector directly onto compatible gas outlet. Push the other end onto the fitting under the nebulizer.
5. Set the gas source to the prescribed flow.
6. Check for secure connections by confirming flow of aerosol through device. If necessary, gently tap the nebulizer until it begins to nebulize. During treatment, periodically tap the nebulizer to minimize residual volume.
7. If provided with a mask, place mask on face with the elastic strap below the ears and around the neck. Gently pull the ends of the strap until the mask is secure. Mold the metal strip on the mask to fit across the bridge of the nose.
8. Shut off gas flow and disconnect supply tubing prior to cleaning. Clean the device after each use.

Disposal instructions:

Discard after 30 days of use or 60 cleaning cycles, whichever comes first. Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Cleaning Guidelines:

These guidelines are for cleaning purposes only and are intended for the removal of medication residue from the nebulizer or any of its parts. Following these guidelines will not sterilize the SVN, mouthpiece, or aerosol mask.

1. Separate the Cap, Jet, and Jar. Do not clean other components.
2. Wash thoroughly with dishwashing soap in warm distilled water.
3. Additionally, you may soak the components in 70% Isopropyl Alcohol for 5 minutes. Rinse the products completely with hot water to remove any residue. Dry thoroughly with a clean cloth or allow to air dry completely before re-assembling.
4. Alternatively, validated hospital cleaning or disinfection procedures may be used.
5. Store in a clean and dry location.

Warnings:

- Do not use with Anaesthetic or Ventilator breathing systems.
- Using a solution other than 0.083% Albuterol Sulfate, in particular a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, as compared to values disclosed below.
- Materials used in the nebulizer may not be compatible with solutions, suspensions, or emulsions that have not been evaluated.
- Non-sterile. Do not sterilize prior to use.
- Do not use if device is soiled or damaged.
- Prior to use, ensure compatibility with gas outlet. Standard connectors are compatible with 5-7mm male fittings and universal connectors are compatible with male DISS 1240 fittings or 5-7mm male fittings.
- Prior to use, ensure all connections are secure. Do not over tighten or force connections. Do not alter device. Improper connections may result in disconnection or leakage during use, which may result in performance degradation and increase the risk of respiratory distress.
- Do not attach the outlet of the nebulizer directly to an endotracheal tube, tracheostomy tube or other device that does not offer an alternate pathway for exhalation.
- Single Patient, Multiple Use Only: Do not reprocess or sterilize. The use of products intended for single patient use on more than one patient creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single patient use only may result in degraded performance or a loss of functionality.
- If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

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CE 2797

Manufactured in Mexico by:
Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093 USA
Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

EC REP Medline International France SAS,
5 rue Charles Lindbergh, 44110 Châteaubriant, France

CH REP **CH rep and importer:** Medline International
Switzerland SARL, c/o MN & Associés SA,
1 Place de Longemalle, 1204 Genève, Switzerland

UK CA 0086 **UK representative and importer in UK:**
Medline Industries Ltd, WA4 6HL, UK

Medline International B.V.,
Nieuwe Stationsstraat 10,
Arnhem, 6811 KS, Netherlands

RA25145 - HNL02195 2025-1