INSTRUCTIONS FOR USE

O.R. Necessities
DYNJLHS1 & DYNJLHS2 Flexible Light Handle Cover

Aspen Surgical Products Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316
Phone 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

EC REP
Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Made in USA

CE
2797

English

Intended use:
Intended to provide a sterile interface between the operating room surgical lamps (non-sterile) and the Surgeon or Nurse (sterile).

Intended User/ Patient Target Groups:
Device is non-patient contacting and is to be used by Healthcare Professionals. Provides a means of manipulating O.R. Lights.

Precaution:
- Verify securement of the light handle cover to the adaptor.
- Do not use if product appears visually damaged or degraded.
- This device is single use only. Reuse of device could result in infection/ contamination and/ or device faliure which could lead to patient harm.
- Do not use product if product is damaged or particulate is found in sterile packaging.

Disposal:
Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Notice to Users and/or Patients in EU:
Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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