0.29 in

0.29 in

HUDSON RCI®

Water Trap

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Instructions for use

Manufactured in Mexico by:
Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093 USA.

CE 2797 EC REP Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

CH REP **CH Rep and importer:** Medline International Switzerland SARL, c/o MN & Associés SA, 1 Place de Longemalle, 1204 Genève, Switzerland

UK CA 0086 **UK representative and importer in UK:** Medline Industries Ltd, WA4 6HL,UK

Medline International B.V., Nieuwe Stationsstraat 10, Arnhem, 6811 KS, Netherlands

Distributed in Australia by:
Medline International Two Australia Pty. Ltd., Sydney, Australia

REF 1650CE

RK24HNL 2024-10
HNL02184

en Water Trap

Intended Use:
Water Traps are intended to trap and drain water that collects in ventilator tubing during respiratory therapy. These devices are intended for use in both acute and non-acute care settings by qualified medical personnel for use with adult populations. Duration of use not to exceed 30 days.

Storage Conditions:
Keep away from sunlight and keep dry. Do not use if the packaging is compromised.

Directions for Use:
- Attach the water trap lid to the jar using a quarter turn clockwise.
- Connect water trap to 22mm tubing. Be sure connections are secure. Place at the lowest point in the circuit to assist in drainage.
- During use if the water does not drain, tap the jar to break the surface tension.
- To empty the trap, turn the jar a quarter turn counterclockwise. The lid will seal the circuit to maintain a closed system. Empty trap when full.

Disposal Instructions:
Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

Warnings:
- Single Use Only. Do not reuse, reprocess, or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
- If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

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