MEDLINE

REF ALARM01

en PATIENT FALLS PREVENTION ALARM MONITOR

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CE UK CA

Manufactured by:
Rondish Company Limited, Unit G&H, 4/F, Block 1, Kwai Tak Ind Ctr 15 - 33 Kwai Tak St Kwai Chung, N.T., Hong Kong

Imported and distributed in Europe by:
Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

RH21ROND 2021-08

en

1. Choose Power Supply
- Remove protective cover from monitor
- Recommended: Use screwdriver to open battery cover and insert 3 x AA alkaline batteries [not included]
- Optional: An AC adapter can be used if necessary – compatible 6V, 300mA [not included]

2. Adjust Monitor Settings
- Adjust the volume to desired level
- Select alarm tone or leave as default
- Set a delay if you see false alarms or only want notification when patient does not return to bed
- The reset button can be disabled using the switch on the left side to prevent patient tampering. When disabled, reset is done using a magnetic fob over the red circle on monitor face

3. Connect Monitor with Sensorpad
- Plug sensorpad cable into connection point at the bottom of the monitor
- Run cable through the strain relief point at bottom of monitor unit

4. Test the Alarm
- Apply pressure on sensorpad – the alarm should sound "beep beep" to indicate it is active
- Status LED will blink every 3 seconds when the system is activated and operating properly
- Release pressure to confirm it alarms
- Press "Reset" button on monitor to silence the alarm

5. Place Sensorpad on Bed or Chair
- Place under shoulder blades for alert when patient sits up in bed
- Place under hips for alert when patient exits the bed or chair

6. Mount the Monitor Unit
- Put soft cover onto the monitor to protect the unit in case of drops
- Attach the monitor to bed or chair frame using buckle strap or wall-mount bracket

TROUBLESHOOTING:

- If you are experiencing false alarms, try setting a 2s delay on the switch inside the battery compartment
- Try disconnecting and reconnecting the sensorpad
- Refer to UDI on side of alarm unit and contact your distributor for assistance

INDICATIONS:

- Store in low humidity (32-85%) warehouse between 10-25 °C
- Battery life: ~6 months in normal usage
- Low battery: Monitor unit will start "chirping" several days before batteries run flat
- Warranty period: 2 years
- This device must be disposed of properly according to local ordinances

Expected safe working lifetime of Monitor: 3 years

WARNING:

Failure to comply with warning may result in injury or death. This device is not suitable for all individuals. Other devices may be required. This device is not a substitute for visual monitoring by caregiver. Be sure cord cannot entangle or choke the user. The manufacturer does not claim that this device will stop elopement and/or stop falls. This device is designed to augment caregivers comprehensive resident mobility management program. Test this device before each use.

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