Suresite®
Transparent Film Dressing with Borders

- de **Transparenter Katheterverband mit Stoffumrandung**
 Gebrauchsanweisung
- fr **Pansement film transparent adhésif à bords renforcés**
 Mode d'emploi
- nl **Transparant filmverband met kleefranden**
 Gebruiksaanwijzing
- it **Medicazione trasparente con bordo**
 Istruzioni per l'uso
- es **Apósito transparente con borde**
 Instrucciones de uso
- pt **Penso de Película Transparente com Moldura**
 Instruções de utilização
- sv **Transparent Filmförband med Ramsystem**
 Bruksanvisning
- fi **Läpinäkyvä kiinnittyvä kalvosidos**
 Käyttöohjeet
- da **Transparent Filmforbinding med Kantsystem**
 Brugsanvisning
- no **Transparent Filmbandasje med Rammesystem**
 Bruksanvisning
- el **Διαφανής κάλυμμα φιλμ με περίγραμμα**
 Οδηγίες χρήσης
- cs **Průhledná fólie s okrajem**
 Návod k použití
- pl **Przezroczysty opatrunek foliowy z ramką**
 Instrukcja stosowania
- sl **Prozorna filmska obloga z robom**
 Navodila za uporabo
- ro **Pansament pe suport de folie transparentă, cu margini**
 Instrucțiuni de utilizare
- bg **Превръзка с прозрачен филм и граници**
 Инструкции за употреба
- et **Äärtega läbipaistev kile**
 Kasutusjuhend
- hu **Átlátszó filmkötszer kerettel**
 Használati utasítás
- lv **Caurspīdīgs plēves pārsteris ar maliņām**
 Lietošanas instrukcija
- lt **Skaidrios plėvelės tvarstis su apvadu**
 Naudojimo instrukcija
- mt **Infaxxatura b'rita trasparenti u bil-borduri**
 Istruzzjonijiet għall-użu
- sk **Priehľadná fólia s okrajmi**
 Návod na použitie
- hr **Prozirna film traka s ivicama za povoj**
 Upute za korištenje
- tr **Kenarlıklı Şeffaf Film Kaplama**
 Kullanım talimatları


Instruction for use

MD STERILE EO

Latex Free, Latexfrei, Sans latex, Latexvrij, Senza lattice, Sin Látex, Sem látex, Latexfri, Lateksiton, Latexfri, Lateksfri, Χωρίς λάτεξ, Bez Latexu, Bezlateksowe, Brez lateksa, Fără latex, Несъдържащ латекс, Lateksivaba, Latex mentes, Nesatur lateksu, Sudėtyje nėra latekso, ma fihx latex, Bez latexu, Bez lateksa, Lateks içermez


REF VAD1610-INT
REF VAD1614-INT
REF VAD1616-INT
REF VAD1633-INT
REF VAD1635-INT
REF VAD1650-INT
REF VAD1655-INT

1


- en Remove first adhesive liner.
- de Entfernen Sie die erste Klebeschicht.
- fr Retirez la première couche d'adhésif.
- nl Verwijder de eerste lijmlaag.
- it Rimuovere il primo liner adesivo.
- es Retire el primer protector de adhesivo.
- pt Remover o primeiro revestimento adesivo.
- sv Avlägsna det första självhäftande skiktet.
- fi Irrota ensimmäisen liimapinnan suojus.
- da Fjern det første klæbende lag.
- no Fjern den første beskyttende plasten.
- el Αφαιρέστε την πρώτη αυτοκόλλητη ταινία.
- cs Odstraňte první lepicí vložku.
- pl Zdjąć samoprzylepną warstwę z jednej strony.
- sl Najprej odstranite lepilno oblogo.
- ro Îndepărtați prima folie adezivă.
- bg Свалете първата адхезивна подложка.
- et Eemaldage esimene kleepuv kate.
- hu Távolítsa el az első ragasztóbetétet.
- lv Noņemiet pirmo adhezīvo aizsargslāni.
- lt Nuplėškite pirmąjį lipnų sluoksnį.
- mt Neħħi l-ewwel kisja li twaħħal.
- sk Odstráňte prvú lepiacu vložku.
- hr Uklonite prvi zaštitni ovoj.
- tr İlk yapışkan astarı çıkarın.

2


- en Apply dressing.
- de Verband anlegen.
- fr Appliquez le pansement.
- nl Breng het verband aan.
- it Applicare la medicazione.
- es Coloque el apósito.
- pt Aplicar o penso.
- sv Applicera förbandet.
- fi Aseta side paikalleen.
- da Påfør bandage.
- no Legg på bandasjen.
- el Εφαρμόστε τον επίδεσμο.
- cs Přiložte krytí.
- pl Nałożyć opatrunek.
- sl Nanesite obvezo.
- ro Aplicați pansamentul.
- bg Поставете превръзката.
- et Kandke side peale.
- hu Helyezze fel a kötést.
- lv Uzlieciet pārsēju.
- lt Užklijuokite tvarstį.
- mt Applika l-garża.
- sk Aplikujte náplasť.
- hr Stavite zavoj.
- tr Sargıyı uygulayın.

3

- en Remove second adhesive liner & smooth edges.
- de Entfernen Sie die zweite Klebefolie und glätten Sie die Kanten.
- fr Retirez la deuxième pellicule adhésive et lissez les bords.
- nl Verwijder de tweede plakfolie en strijk de randen glad.
- it Rimuovere il secondo liner adesivo e lisciare i bordi.
- es Retire el segundo protector de adhesivo y alise los bordes.
- pt Remover o segundo revestimento adesivo e as bordas lisas.
- sv Avlägsna det andra självhäftande skiktet och jämna till kanterna.
- fi Irrota toisen liimapinnan suojus ja silota reunat.
- da Fjern det andet klæbende lag og udglat kanterne.
- no Fjern den andre beskyttende plasten og jevn ut kantene.
- el Βγάλτε τη δεύτερη αυτοκόλλητη ταινία και λειάνετε τα άκρα.
- cs Odstraňte druhou lepicí vložku a uhlaďte okraje.
- pl Usunąć samoprzylepną warstwę z drugiej strony i wygładzić krawędzie.
- sl Odstranite drugo lepilno oblogo in zgladite robove.
- ro Îndepărtați a doua folie adezivă și neteziți marginile.
- bg Свалете втората адхезивна подложка и изгладете ръбовете.
- et Eemaldage teine kleepuv kate ja siluge servad.
- hu Távolítsa el a második ragasztóbetétet és simítsa el a széleket.
- lv Noņemiet otro adhezīvo aizsargslāni un nogludiniet malas.
- lt Nuplėškite antrąjį lipnų sluoksnį ir išlyginkite kraštus.
- mt Neħħi t-tieni kisja li twaħħal u t-truf lixxi.
- sk Odstráňte druhú lepiacu vložku a vyhlaďte okraje.
- hr Skinite drugi zaštitni ovoj i zagladite rubove.
- tr İkinci yapışkan astarı çıkarın ve kenarları yumuşatın.

4

- en Apply securement tape strip.
- de Bringen Sie den Fixierstreifen an.
- fr Appliquez la bande de sécurisation.
- nl Breng bevestigingstape strip aan.
- it Applicare la striscia di nastro di fissaggio.
- es Coloque la cinta de fijación.
- pt Aplicar fita adesiva de segurança.
- sv Applicera fasthållningstejp.
- fi Kiinnitä teippi.
- da Påfør sikringstapestrimmel.
- no Fest med teipstrimmelen.
- el Εφαρμόστε ταινία ασφαλείας.
- cs Aplikujte upevňovací pásku.
- pl Zastosować pasek taśmy zabezpieczającej.
- sl Nanesite lepilni trak za pritrjevanje.
- ro Aplicați banda adezivă de siguranță.
- bg Поставете лентата за сигурно фиксиране.
- et Pange peale kinnitusriba.
- hu Helyezze fel a rögzítő szalagcsíkot.
- lv Uzlieciet nostiprinošās līmlentes loksni.
- lt Užklijuokite apsauginę juostą.
- mt Applika strixxa tejp biex iżommha f'postha.
- sk Použite fixačnú pásku.
- hr Stavite traku za pričvršćivanje.
- tr Sabitleyici bant askıyı uygulayın.

en

Indications:
Intended to secure and maintain catheters. Can be left in place for up to 7 days.

General Precautions & Observations:
Can be used as a primary (stand alone) or secondary dressing (with one fixation method) .
Do not stretch the film during application.
Should any signs of redness or irritation appear, discontinue use and consult a healthcare professional.

Contraindications:
· Do not use on damaged skins or wounds
· Those devices are not wound dressing
· Must not be used on neonates or burn patients

Warning:
Using the Suresite Transparent Film Dressing with Borders to cover a drug delivery patch could lead to faster absorption, increasing the risk of overdose.
Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient.
If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Disposal instructions:
Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

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Medline International France SAS,
UK CA 0086 **UK representative and importer in UK:** Medline Industries Ltd, WA4 6HL, UK