
CENTURION®
Medical Products

Compass® Lumbar Puncture

English....2
German....6
French....10
Dutch....14
Italian....18
Spanish....22
Portuguese....26
Swedish....30
Czech....34
Slovak....38
Slovenian....42
Norwegian....46
Polish....50
Finnish....54
Danish....58
Greek....62
Romanian....66
Bulgarian....70
Estonian....74
Hungarian....78
Latvian....82
Lithuanian....86
Croatian....90
Russian....94
Turkish....98
Korean....102
Arabic....106
Hebrew....110
Malay....114
Serbian....118

English

CENTURION®
Medical Products

Compass® Lumbar Puncture

Intended Use
The Compass® Lumbar Puncture disposable pressure transducer with integrated digital display is intended for direct measurement and monitoring of intracranial pressure. Compass devices are intended for use by a licensed healthcare professional within the general population including pediatrics.

Indications for Use
The Compass® Lumbar Puncture is indicated for use in the measurement and monitoring of intracranial pressure.

Caution: Federal (USA) law restricts this device to sale by or on on the order of a physician.

Information for Use: Please review prior to use. Sterilized using ethylene oxide. Device is sterile if pouch is unopened and undamaged. This product is designed for single use only. Do not resterilize or reuse. Not for use in oxygen-rich environments.

WARNING: No modification of this equipment is allowed. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

The Compass® device contains a Lithium Manganese button cell battery that cannot be removed from the device housing. The whole Compass® device should be disposed of as medical waste as facilitated by the healthcare facility's protocol and local member state and/or government requirements.

A Suggested Procedure

1. Remove Compass® from its sterile package.
2. Activate Compass® by depressing power button until display reads "00" (refer to Fig. 1).

WARNING: Turn on Compass® prior to attaching needle to avoid incorrect baseline pressure.

3. Remove cap (refer to Fig. 2).
4. Once the spinal needle and stylet are inserted through the skin and subcutaneous tissue[1], remove stylet and attach Compass® to needle (refer to Fig 3). If needed (e.g. during procedures with patient in prone position), extension tubing can be used to position the Compass® at the desired location (level of the spine) to obtain accurate pressure reading.

WARNING: To ensure physiologic pressure measurements, completely cover rear Compass® port with thumb. Thumb may be removed momentarily to zero device.

5. Cover rear Compass® port with thumb and advance spinal needle into CSF space per hospital protocol; pressure should increase upon entry of spinal needle into CSF space.
6. Note opening pressure (refer to Fig. 4).


WARNING: Once needle has entered CSF space, wait until pressure has stabilized (2-3 seconds) before recording opening pressure.

7. Collect CSF sample from proximal port of Compass® (refer to Fig. 5).


 Alternate Technique: Compass® can be removed from spinal needle and CSF sample collected from spinal needle hub.

WARNING: Do not reinfuse CSF or other bodily fluid into the patient after any such fluid is drawn into the Compass®.


1


2


3

4

5

WARNING: Do not reuse the Compass® for more than one procedure in the same patient. Use a new Compass® device for each procedure. Reuse may lead to patient infection or cross contamination due to lack of sterility.

8. Remove Compass® from spinal needle, replace stylet, and withdraw spinal needle[2].

References:


1. Nigrovic, L.E., Kuppermann, N. & Neuman, M.I. Risk factors for traumatic or unsuccessful lumbar punctures in children. Ann Emerg Med 49, 762–71 (2007).

2. Williams, J., Lye, D.C. & Umpathi, T. Diagnostic lumbar puncture: minimizing complications. Intern Med J 38, 587–91 (2008).

Device Power Off

Press button three times in a ("PUSH","PUSH","HOLD") sequence; quickly push button twice and then hold third button push until device turns off. Turn device back on to reset.

Note: If Compass® pressure reading does not return to zero after removal from tissue, turn off device to reset the zero.

Compass® Lumbar Puncture Specifications

Characteristics	Value
Pressure Range	-199 to +199 cm H_2O
Accuracy	±2.7 cm H_2O below 68 cm H_2O ±3% of reading above 68 cm H_2O
Zero Drift	±2 cm H_2O per 4 hours
Operating Life	8 hours
Operating Temperature	15 to 40°C (59 to 104°F)
Storage Temperature	-20 to 49°C (-4 to 120°F)
Light Sensitivity	< 2 cm H_2O
Lithium Battery Voltage	3 Volts

Software Error Codes

If the software detects device malfunction, the following codes will be displayed and the device rendered inoperable.

Code	Meaning
SF	Software Failure
HF	Hardware Failure
PF	Pressure Sensor Failure
bF	Battery Failure-Less than 30 minutes of use remains at device power – on
EE	Environmental Error-Outside operating temperature range
CE	Calibration Error-Pressure values outside required accuracy (±3 cm H_2O)

©2024 Medline Centurion
Centurion Medical Products, LP
100 Centurion Way
Williamston, MI 48895 USA
www.medline.com
1-800-MEDLINE

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EC REP Medline International France SAS
Medline International B.V.
BSI Assurance UK Ltd,
Medline Industries Ltd,