MEDLINE

HUDSON RCI®

Tracheostomy Adaptor

de **Tracheotomie-Adapter** - Gebrauchsanweisung
fr **Adaptateur pour trachéotomie** - Mode d'emploi
nl **Tracheotomie-adapter** - Gebruiksaanwijzing
it **Adattatore per tracheostomia** - Istruzioni per l'uso
es **Adaptador para Traqueostomía** - Instrucciones de uso
pt **Adaptador de Traqueostomia** - Instruções para Uso
sv **Trakeostomiadapter** - Instruktioner för användning
fi **Trakeostomiasovitin** - Käyttöohjeet
da **Trakeostomiadapter** - Instruktion
no **Trakeostomiadapter** - Bruksanvisning
el **Προσαρμογέας Τραχειοτομής** - Οδηγίες χρήσης
cs **Tracheostomický adaptér** - Návod k použití
pl **Adapter do tracheostomii** - Instrukcja obsługi
sl **Adapter za traheostomijo** - Navodila za uporabo
ro **Adaptor de traheostomie** - Instrucțiuni de folosire
bg **Адаптер за трахеостома** - Инструкции за употреба
hu **Trakeostomiás adapter** - Használati utasítás
sk **Tracheostomický adaptér** - Návod na použitie
tr **Trakeostomi Adaptörü** - Kullanim talimatlari

MD

Mandatory action to read instructions for use. Unbedingt Gebrauchsanweisung beachten, Obligation de lire la notice d'utilisation, Verplichte actie om de gebruiksaanwijzing te lezen, Azione obbligatoria: leggere le istruzioni per l'uso, Acción obligatoria de leer las instrucciones de uso, Acção obrigatória para ler as instruções de utilização, Det är obligatoriskt att läsa bruksanvisningen, Käyttöohjeiden lukeminen on pakollista, Det er påkrævet at læse brugsanvisningerne, Obligatorisk tiltak som må leses før bruk, **Η συσκευασία των οδηγιών είναι υποχρεωτική πριν από τη χρήση**, Přečtěte si návod k použití, Należy koniecznie przeczytać instrukcję obsługi, Obvezni ukrep je branje navodil za uporabo, Citirea instrucțiunilor de utilizare este obligatorie, Задължително е да се прочетат инструкциите за употреба, Kötelező intézkedés a használati utasítás elolvasása, Povinné přečítanie si návodu na použitie, Kullanım talimatlarını okumak zorunludur

Latex Free, Latexfrei, Sans latex, Latexvrij, Senza lattice, Sin Látex, Sem látex, Latexfri, Lateksiton, Latexfri, Latekszfri, Χωρίς λατέξ, Bez Latexu, Bezlateksowe, Brez lateksa, Fără latex, Несъдържащ латекс, Latex mentes, Bez latexu, Latéks ücretsiz

Instruction for use

REF 241501CE

en Instructions for Use Tracheostomy Adaptor

Intended Use:
Tracheostomy Adaptor is used in conjunction with the tracheostomies in order to provide a comfortable secure connection for oxygen therapy. It allow for the safe and effective delivery of oxygen therapy to the patients trachea.

Indication for Use:
It is indicated to be used by spontaneously breathing tracheostomized patients only.

Product Configuration:
Tracheostomy Adaptor consists part as below:
I. 15mm (Female) connector with Cap Deflector
II. Lanyard
III. Swivel Connector (15F/22M) with tubing

Connection
One (1) 15mm I.D., one (1) 15mm I.D./22mm O.D., comply to ISO 5356-1.
Note: Device to be used only under supervision of a physician.

Direction of Use
1. Connect 15mm ISO connection to tracheostomy tube.
2. Connect 15/22mm ISO connection to oxygen therapy delivery.
3. Adjust the lanyard accordingly for maximum comfort andsecurement. Optional: Attach the lanyard clip to the patient's gown or clothing for better securement.
4. Ensure all connections are secure.

Warning:
- Do not obstruct exhalation ports.
- Check for condensate regularly. Drain as required.
- Do not use with an air entrainment device.
- Continuously monitor patient condition to ensure adequate flow of gas is being delivered to patient.
- Do not stretch the tubing to remove condensation, damage or malfunction may occur.
- Do not use if package is broken or product is damaged.
- Ensure that all connections are secure and safe.
- To ensure loading on tracheostomy tube is kept to a minimum, make sure lanyard is fitted appropriately.
- Use of lanyard may not be suitable for patients with damaged skin.
- The unit should be replaced immediately if there is increasing resistance due to accumulation of secretions.

Cautions
Warning: Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient. If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

Possible Complication
The device might result in higher flow resistance due to accumulation of secretions.

Disposal
Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologicallyhazardous waste.

Storage Instruction
Keep away from sunlight and keep dry. Do not use if the packaging is compromised.

Useful Life
30 This product is intended for use for a maximum of 30 days.

Explanation of important symbols and markings on the product label

Qty	Quantity, Number of items
	Not made with natural rubber latex
	Do not use if package is damaged

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