HUDSON RCI®

Tracheostomy Adaptor

REF 2415-01


CE 0124

Teleflex Medical Sdn. Bhd.
Lot PT 2577, Jalan Perusahaan 4,
34600 Kamunting, Perak, Malaysia
Made in Malaysia

© Copyright 2017 Teleflex

EC REP

Teleflex Medical
IDA Business & Technology Park,
Dublin Road, Athlone,
Co. Westmeath, Ireland

en

Instructions for Use
Tracheostomy Adaptor

Intended Use:
Tracheostomy Adaptor is used in conjunction with the tracheostomies in order to provide a comfortable secure connection for oxygen therapy. It allow for the safe and effective delivery of oxygen therapy to the patients trachea.

Indication For Use:
It is indicated to be used by spontaneously breathing tracheostomized patients only.

Product Configuration:
Tracheostomy Adaptor consists part as below:
I. 15mm (Female) connector with Cap Deflector
II. Lanyard
III. Swivel Connector (15F/22M) with tubing

Connection
One (1) 15mm I.D., one (1) 15mm I.D./22mm O.D., comply to ISO 5356-1.

Note: Device to be used only under supervision of a physician.

Direction of use
1. Connect 15mm ISO connection to tracheostomy tube.
2. Connect 15/22mm ISO connection to oxygen therapy delivery.
3. Adjust the lanyard accordingly for maximum comfort and securement. Optional: Attach the lanyard clip to the patient's gown or clothing for better securement.
4. Ensure all connections are secure.

⚠ Warning:
- Do not obstruct exhalation ports.
- Check for condensate regularly. Drain as required.
- Do not use with an air entrainment device.
- Continuously monitor patient condition to ensure adequate flow of gas is being delivered to patient.
- Do not stretch the tubing to remove condensation, damage or malfunction may occur.
- Do not use if package is broken or product is damaged.
- Ensure that all connections are secure and safe.
- To ensure loading on tracheostomy tube is kept to a minimum, make sure lanyard is fitted appropriately.
- Use of lanyard may not be suitable for patients with damaged skin.
- The unit should be replaced immediately if there is increasing resistance due to accumulation of secretions.

Cautions

Warning: Single use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

DEHP

Contains DEHP
The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical devices containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Possible Complication
The device might result in higher flow resistance due to accumulation of secretions.

Disposal
The product must be disposed of in accordance with hospital protocol /state regulations.

Storage Instruction
Keep away from sunlight and keep dry. Do not use if the packaging is compromised.

Useful Life

30 This product is intended for use for a maximum of 30 days.

Explanation of important symbols and markings on the product label.

Qty	Quantity, Number of items
LATEX	Not made with natural rubber latex
	Do not use if package is damaged

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