



OPEN HERE



STERILE URINARY CATHETER FOR SINGLE USE (Foley catheter)
CATETERE URINARIO STERILE MONOUSO

WARNING: Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

ATTENZIONE: Non utilizzare unguenti o lubrificanti a base di petrolio. Danneggeranno il catetere e faranno scoppiare il palloncino.

DIRECTIONS: After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

INDICAZIONE: Dopo il cateterismo asettico gonfiare il palloncino di ritenzione inserendo la punta della siringa nella valvola e introdurre la quantità necessaria di acqua sterile.



Jiangsu Suyun Medical Materials Co., Ltd.
No. 18 Jin Qiao Road Qiao Industrial Park 222002 Lianyungang, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
Tel: +86-518-85503111 Fax: +86-518-85460033 Made in China

tel.	nr.	Shanghai International Fidelity Corp. GmbH (Europe) Eiffelsaal 83, 20087 Hamburg, GERMANY
------	-----	--

 **Gima S.p.A.**
Via Marconi, 1
20060 Gessate (MI) - Italy
www.gimaitaly.com

978-1-107-01033-0

Type: 2-WAY-STANDARD TIP
Size: 4.0mm(12Fr/CH). 30cc/ml

REF S5122.12 (GIMA 22064)

LOT 000014



240403

290401

☐ [Download](#)
☐ [View](#)
☐ [Print](#)
☐ [Share](#)